

Spor-Klenz[®] Concentrate

Sporicidal Agent



Spor-Klenz Concentrate is a liquid, cold sterilant/disinfectant specifically formulated for sterilization and disinfection of hard surfaces. Spor-Klenz Concentrate is a stabilized blend of peracetic acid, hydrogen peroxide and acetic acid, which provides fast, effective microbial control, including spores. Spor-Klenz Concentrate offers:

- *Microbial control against bacteria, fungi, viruses and bacterial spores*
- *Versatility of use as a sterilant, high-level disinfectant or sanitizer*
- *0.22 micron filtration*

Applications

- Aseptic areas
- Critical Pharmaceutical, Biotechnology, and Medical Device Manufacturing areas
- Processing areas including:
 - Equipment
 - Hoods
 - Countertops
 - Laminar Flow Curtains
- Other FDA Regulated facilities

Dilution Rate/Ordering Information

Dilution Rate:	1:99
1-quart bottle x 4 per case:	6520-26

Features and Benefits

- Filtered to 0.22 micron to reduce particulate levels
- Each quart makes 25 gallons of use dilution, a cost-effective method
- Easy, one-step dilution eliminates the need for multiple products to deliver a sterilant/high-level disinfectant
- Use dilution has a shelf life of up to seven days, which reduces product waste due to expiration
- Low Toxicity Profile to support worker safety
- Certificate of Analysis, Filtration and Sterility with each case to support pharmaceutical documentation requirements
- Complete lot traceability
- Sterilant and high-level disinfectant uses which provide flexible options depending on facility needs

- Sterilant	11 hrs.	@ 20° C
- Broad Spectrum Disinfectant	10 min.	@ 20° C



1. PRODUCTS AT THE BASE OF THE PYRAMID ARE MOST FREQUENTLY USED AND ARE GENERALLY NOT SPORICIDAL.
2. PROGRESSION UP THE PYRAMID INDICATES STRONGER PERFORMANCE OVERALL AND A BROADER SPECTRUM OF CLAIMS.

STERIS Corporation is a global provider of contamination control systems, products and services for the pharmaceutical, biotech, medical device, laboratory research and industrial markets. Our solutions cover all aspects of product discovery and development – across research lab, animal research lab and production/manufacturing processing.

STERIS Offices Worldwide

Benelux	32 2 523 2488
Canada	800 661 3937
France	33 1 44 882688
Germany	49 2233 69990
Hong Kong	852 2 563 3623
Italy	39 02 66 80 53 10
Japan	81 78 321 2271
Latin America	1 305 442 8202
Nordic	358 9 25851
Singapore	65 841 7677
Spain	34 916 585 920
Sweden	46 152 228 30
United Kingdom	44 1256 840400

Publication ID #M1473EN.0107, Rev. B
 ©2001 by STERIS Corporation.
 All rights reserved.

Technologies to Prevent Infection and Contamination



STERIS Corporation
 Scientific Division
 5960 Heisley Road
 Mentor, OH 44060-1834 ■ USA
 440-354-2600 ■ 800-444-9009
 www.steris.com