


The Vesphene III Disinfectants are specifically formulated for the routine disinfection of hard, non-porous surfaces. The Vesphene III Disinfectants are processed and packaged differently for specific applications. At a use dilution of rate of 0.8% v/v (1:128), the Vesphene III Disinfectants kill a broad spectrum of common bacteria and fungi. The alkaline Vesphene III Disinfectant formulation can be used alone or in rotation with the LpH® III Disinfectants. These products have an extensive documentation package to meet your validation objectives and are supported by the STERIS Technical Services team for individualized application and validation assistance.



Vesphene® III Disinfectants

Our latest high pH phenolic disinfectant

APPLICATIONS

- Surface disinfection and microbial control of:
 - Floors
 - Walls
 - Ceilings
 - Countertops
 - BSL hoods
 - Carts and other equipment
- Products may be applied by:
 - Mopping
 - Wiping
 - Spraying

FEATURES

Proven tuberculosis (*Mycobacterium tuberculosis*, var. bovis) efficacy

Broad spectrum efficacy against bacteria, fungi, and viruses using AOAC protocols

Superior material compatibility

Formulated with biodegradable surfactants

Active ingredients are ortho phenylphenol (OPP) and ortho benzyl para chloro phenol (OBPCP)

Sterile version is gamma irradiated and each lot USP sterility tested

Manufactured in accordance with cGMP with complete lot traceability

BENEFITS

Indicative of potential efficacy against less challenging microorganisms

Confirms efficacy against many common microbes


Allows for high frequency use on most surfaces

Reduces environmental impact of disposal while maintaining excellent cleaning ability

Registered by the United States Environmental Protection Agency (EPA)

Introduces minimal bioburden into your facility

Supports quality assurance and quality control



The foundation of any contamination control program is the disinfectant that should be used daily. Disinfectants should be used on a scheduled rotation with a sporicide based on your environmental data. Historically, multiple disinfectants have been rotated, but if your disinfectant has broad spectrum efficacy, this is not required. Either used alone or rotated with the acidic LpH III Disinfectant family, the alkaline Vesphene III Disinfectants offer excellent material compatibility for daily use. Our Technical Service group can help customize your contamination control program to take advantage of this excellent product.

¹Products that fall into the categories at the bottom of the pyramid are most frequently used and are generally not sporicidal. Progression up the pyramid indicates stronger performance overall and a broader spectrum of claims.

ORDERING INFORMATION

Vesphene III se Disinfectant

Filtered for consistent quality for use in a variety of areas

United States

1 gallon bottle	4 per case	6475-08
5 gallon pail	Each	6475-05
55 gallon drum	Each	6475-01

Canada

3.78 L bottle	4 per case	6475-08CA
18.93 L pail	Each	6475-05CA
208.18 L drum	Each	6475-01CA

Vesphene III st Disinfectant

Double bagged and gamma irradiated through a validated process for easy introduction into your aseptic areas

United States

2 fluid ounce pouch	144 per case	1S20-29
1 gallon bottle	4 per case	1S20-08

Canada

59.1 mL pouch	144 per case	1S20-29CA
---------------	--------------	-----------

One Solution st® Vesphene III Premeasured WFI

A pouch of disinfectant attached to a gallon of USP Water For Injection and double bagged and gamma irradiated through a validated process for easy introduction into your aseptic areas

United States

1 fluid ounce pouch and 1 gallon of USP WFI	4 per case	1S15-CX
---	------------	---------



STERIS CORPORATION
5960 Heisley Road
Mentor, OH 44060, USA
T +1 800 444 9009

www.STERISLifeSciences.com

STERIS Offices Worldwide

Asia Pacific	+65 (0) 68 41 7677
Canada	+1 800 444 9009
Finland	+358 9 25851
France	+33 (0) 2 38 84 85 40
Germany	+49 (0) 221 466120
Italy	+39 (0) 2 21303 424
Latin America	+1 800 444 9009
Spain	+34 (0) 916 585 920
United Kingdom	+44 (0) 1256 840400

STERIS Life Sciences

For over 100 years, STERIS has been known as a global leader, trusted partner and solutions provider in the field of sterilization and contamination control. Today, STERIS continues building on this heritage by advancing the Science of sterilization, cleaning and infection control while offering Solutions that meet our Customers' needs and high standards. STERIS is dedicated to helping you enhance the Life of your patients and the life of your equipment. From Formulated Chemistries, to Capital Equipment, and to Parts and Services, STERIS Life Sciences is your Solutions provider.

© 2017 STERIS Corporation
All rights reserved

No part of this publication may be reproduced in any material form without the permission of STERIS Corp.